

# UGT1A10 Polyclonal Antibody



Catalog Number: 166609

## Product name

UGT1A10 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to UGT1A10

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human UGT1A10 (NP\_061948.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - UGT1A10 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using UGT1A10 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.