

# STMN3 antibody



Catalog Number: 107609

## Product name

STMN3 antibody

## Specificity

Human; other species not tested.

## Antibody description

STMN3 Mouse Monoclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of STMN3 recombinant protein (Accession Number: NM\_015894). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107609(STMN3 antibody) at dilution of 1:1000