

Anti-FH/Fumarase Rabbit antibody



Catalog Number: 168609

Product name

Anti-FH/Fumarase Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

FH/Fumarase Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human FH

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/500

ICC/IF: 1/200

FC: 1/100

IP: 1/20

Validations



Western blot detection of FH/Fumarase in Jurkat,Rat Brain,C6,CHO-K1 cell lysates using FH/Fumarase Rabbit mAb(1:1000 diluted).Predicted band size:55kDa.Observed band size:49kDa.