

EED antibody

Catalog Number: 110309

Product name

EED antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EED Rabbit Polyclonal antibody. Positive IP detected in PC-3 cells. Positive WB detected in COLO 320 cells, HeLa cells, K-562 cells, PC-3 cells. Observed molecular weight by Western-blot: 50-68 kDa

Preparation

This antibody was obtained by immunization of EED recombinant protein (Accession Number: NM_001308007). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110309(EED antibody) at dilution of 1:1000



IP Result of anti-EED (IP:Catalog No:110309, 4ug; Detection:Catalog No:110309 1:500) with PC-3 cells lysate 3200ug.