

PTN antibody

Catalog Number: 114309

Product name

PTN antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PTN Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 15-18 kDa

Preparation

This antibody was obtained by immunization of PTN recombinant protein (Accession Number: NM_002825). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114309(PTN antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:114309(PTN Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:114309(PTN Antibody) at dilution of 1:50

PTN antibody



Catalog Number: 114309

(under 40x lens)