

# Rabbit anti HA-Tag pAb



Catalog Number: 166309

## Product name

Rabbit anti HA-Tag pAb

## Antibody description

Polyclonal antibody to Rabbit anti HA-Tag

## Preparation

Antigen: A synthetic peptide of YPYDVPDYA (HA-tag)-KLH

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF, ChIP


## Dilutions

WB 1:500 - 1:1000

IF 1:50 - 1:100

ChIP 1:20 - 1:50

## Validations



Western blot - Rabbit anti HA-Tag pAb

Western blot analysis of over-expressed HA-tagged protein in 293T cell using HA-tag antibody at different dilution. Each lane was loaded with 2 ug cell lysate. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Blocking buffer: 3% nonfat dry milk in TBST.