

# Anti-SLC5A2 antibody



Catalog Number: 176309

## Product name

Polyclonal

Anti-SLC5A2 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to SLC5A2

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human slc5a2. aa range:101-150

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1:: Western blot analysis of the lysates from Jurkat cells using SLC5A2 antibody.