

Anti-SLC5A2 antibody



Catalog Number: 176309

Product name

Anti-SLC5A2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to SLC5A2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human slc5a2. aa range:101-150

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of the lysates from Jurkat cells using SLC5A2 antibody.