

GST tag antibody



Catalog Number: 107309

Product name

GST tag antibody

Specificity

Tag

Antibody description

GST tag Mouse Monoclonal antibody. Positive WB detected in recombinant protein. Positive IP detected in Recombinant protein protein.

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: FN315687). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:5000-1:50000

IP: 1:5000-1:50000

Validations



Western blot of GST recombinant protein with anti-GST-tag (Catalog No:117321) at various



dilutions.

IP Result of anti-GST Tag (IP:Catalog No:117321, 4ug; Detection:Catalog No:117321 1:20000) with Recombinant protein protein lysate 1000ug.