

Mouse Activin A / INHBA (His Tag) recombinant protein



Catalog Number: 500709

General Information

Gene Name Synonym

Activin beta-A chain

Protein Construction

A DNA sequence encoding the pro form of mouse INHBA (NP_032406.1) (Met 1-Ser 424) was expressed, with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

1. Measured by its binding ability in a functional ELISA.
2. Immobilized human ACVR2B (Cat:503817) at 10 µg/mL (100 µl/well) can bind biotinylated mouse INHBA-His (Cat:500709), The EC₅₀ of biotinylated mouse INHBA-His (Cat:500709) is 0.161 µg/mL.
3. Immobilized mouse INHBA-his at 10 µg/mL (100 µl/well) can bind human Follistatin Protein, The EC₅₀ of human Follistatin Protein is 0.39 µg/mL.

Purity

> 88 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 21

Molecular Mass

The secreted recombinant mouse INHBA comprises 415 amino acids and has a calculated molecular mass of 46.5 kDa. As a result of glycosylation, the apparent molecular mass of rmINHBA is approximately 18 & 43 kDa corresponding to the mature inhibin beta A chain and the pro domain respectively in SDS-PAGE under reducing conditions.

Formulation

- Lyophilized from sterile PBS, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
 2. Please contact us for any concerns or special requirements.

Usage Guide

Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution


Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE

Mouse Activin A / INHBA (His Tag) recombinant protein



Catalog Number: 500709



Mouse Inhibin beta A / INHBA / Activin A Protein (His Tag) SDS-PAGE