

KCTD14 antibody



Catalog Number: 112019

Product name

KCTD14 antibody

Specificity

human, mouse; other species not tested.

Antibody description

KCTD14 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse brain tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of KCTD14 recombinant protein (Accession Number: NM_001282406). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112019(KCTD14 antibody) at dilution of 1:400