

PM20D2 antibody

Catalog Number: 114019

Product name

PM20D2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PM20D2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IHC detected in human placenta. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of PM20D2 recombinant protein (Accession Number: NM_001010853). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114019(PM20D2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human placenta using Catalog No:114019(PM20D2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:114019(PM20D2 antibody) at dilution of 1:50 (under 40x lens)