

Anti-FAT4 antibody



Catalog Number: 175019

Product name

Anti-FAT4 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to FAT4

Preparation

This antigen of this antibody was synthetic peptide within human fat tumor suppressor homolog 4 aa 2260-2313.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC

Dilutions

WB:1:500-1:1,000

ICC: 1:50-1:100

Validations




Fig1: Western blot analysis on recombinant FAT4 using anti-FAT4 polyclonal antibody.




Fig2: ICC image of FAT4 antibody stained D3 (Murine Embryonic Stem Cell Line).The secondary antibody (red) was goat anti-rabbit IgG (H+L) alexa555 conjugated.