

# GSTM1 Polyclonal Antibody



Catalog Number: 166019

## Product name

GSTM1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to GSTM1

## Preparation

Antigen: A synthetic peptide of human GSTM1

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, FC


## Dilutions

WB 1:200 - 1:500

IHC 1:50 - 1:100

FC 1:20 - 1:50

## Validations



Western blot - GSTM1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GSTM1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.