

Anti-IGFBP4 antibody



Catalog Number: 106019

Product name

Anti-IGFBP4 antibody

Immunogen

[Cynomolgus IGFBP4 \(His Tag\) recombinant protein](#)

Specificity

Cynomolgus IGFBP4

Antibody description

Rabbit polyclonal to IGFBP4

Preparation

Produced in rabbits immunized with purified, recombinant Cynomolgus IGFBP4 (rh IGFBP4; ; Met1-Glu258). IGFBP4 specific IgG was purified by Cynomolgus IGFBP4 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Cynomolgus IGFBP4. The detection limit for Cynomolgus IGFBP4 is < 0.039 ng/well.


IHC-P: 0.1-2 µg/mL

Validations



IGFBP4 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of cynomolgus IGFBP4 in cynomolgus stomach with rabbit polyclonal antibody (0.3 µg/mL, formalin-fixed paraffin embedded sections).



IGFBP4 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Anti-IGFBP4 antibody



Catalog Number: 106019

Immunochemical staining of cynomolgus IGFBP4 in cynomolgus kidney with rabbit polyclonal

antibody (0.3 µg/mL, formalin-fixed paraffin embedded sections).