

# Anti-GSC2 antibody



Catalog Number: 176019

## Product name

Anti-GSC2 antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to GSC2

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gsc2. aa range:131-180

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations



Fig1;; Western Blot analysis of various cells using GSC2 Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).