

# EXTL3 antibody

Catalog Number: 110419

## Product name

EXTL3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EXTL3 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of EXTL3 recombinant protein (Accession Number: XM\_047421518). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110419(EXTL3 antibody) at dilution of 1:100



(under 25x lens)

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110419(EXTL3 Antibody) at dilution of 1:600