

# MCM2 antibody



Catalog Number: 107419

## Product name

MCM2 antibody

## Immunogen

[Human MCM2 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

MCM2 Mouse Monoclonal antibody. Positive WB detected in K-562 cells, HEK-293 cells. Observed molecular weight by Western-blot: 120-125 kDa

## Preparation

This antibody was obtained by immunization of MCM2 recombinant protein (Accession Number: NM\_004526). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107419(MCM2 Antibody) at dilution of 1:2000