

# Anti-Hexokinase I (7A7) Mouse antibody



Catalog Number: 168419

## Product name

Anti-Hexokinase I (7A7) Mouse antibody

## Specificity

Human, Rat, Mouse

## Antibody description

Hexokinase I (7A7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human HK1 expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P, ICC/IF

## Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

ICC: 1/200 - 1/1000

ELISA: 1/10000

## Validations



Western blot detection of Hexokinase 1 in Jurkat, HepG2 and HeLa cell lysates using Hexokinase 1 mouse mAb (dilution 1:500). Predicted band size: 102KDa. Observed band size: 102KDa.