

Human IFNGR1 (His tag) recombinant protein



Catalog Number: 518419

General Information

Gene Name Synonym

Interferon gamma receptor
1,IFNGR1,CDw119,CD119

Protein Construction

A DNA sequence encoding the Human IFNGR1
Glu18-Ser245 is expressed with a 6His tag at the
C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM
PB,150mM NaCl,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human IFNGR1 (His tag) recombinant protein