

# c-Myc antibody



Catalog Number: 109419

## Product name

c-Myc antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

c-Myc Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HeLa cells, NCCIT cell. Observed molecular weight by Western-blot: 62-65 kDa, 50 kDa

## Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM\_002467). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:109419(MYC antibody) at dilution of 1:800