

# Anti-KCNN2 antibody

Catalog Number: 175819

## Product name

Anti-KCNN2 antibody

## Specificity

Human, Rat, Mouse

## Antibody description

Rabbit polyclonal antibody to KCNN2

## Preparation

This antigen of this antibody was synthetic peptide of kcnn2 (sk2)

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:1000-1:2000

IHC-P: 1:100-200

## Validations




Fig1:: Western blot analysis of 1) Rat BrainTissue, 2) Mouse Brain Tissue, 3) HepG2 with KCNN2(SK2) Rabbit pAb diluted at 1:2,000.




Fig2:: Immunohistochemical analysis of paraffin-embedded Human BrainTissue using KCNN2(SK2) Rabbit pAb diluted at 1:200.