

# PPEF1 antibody



Catalog Number: 114119

---

## Product name

PPEF1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PPEF1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, human brain tissue. Observed molecular weight by Western-blot: 76 kDa

## Preparation

This antibody was obtained by immunization of PPEF1 recombinant protein (Accession Number: XM\_005274553). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114119(PPEF1 antibody) at dilution of 1:300