

Anti-APOC3 antibody

Catalog Number: 175119

Product name

Anti-APOC3 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to APOC3

Preparation

This antigen of this antibody was recombinant protein with human apoc3 aa 1-100.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol. Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC

Dilutions

WB: 1:500-1:2000

IHC-P: 1:50-1:200

FC: 1:50-1:100

Validations




Fig1: Western blot analysis of APOC3 on human liver tissue lysate using anti-APOC3 antibody at 1/500 dilution.




Fig2: Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-APOC3 antibody. Counter stained with hematoxylin.




Fig3: Flow cytometric analysis of HepG2 cells with APOC3 antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti-rabbit IgG was used as the secondary antibody.