

# Anti-STAT3 (4C8) Mouse antibody



Catalog Number: 168119

## Product name

Anti-STAT3 (4C8) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human, Mouse, Rat, Monkey, Hamster

## Applications

WB

## Antibody description

STAT3 (4C8) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human STAT3 protein fragments expressed in E.coli

## Validations

## Formulation


Ascites containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of STAT3 in CHO-K1, C6, 3T3, K562, CEM, Jurkat, HeLa, MCF7, COS7, 293T, A431 and A549 cell lysates using STAT3 mouse mAb (1:1000 diluted). Predicted band size: 88KDa. Observed band size: 88KDa.