

# HOMER3 antibody



Catalog Number: 111519

---

## Product name

HOMER3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HOMER3 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 40kd, 45kd

## Preparation

This antibody was obtained by immunization of HOMER3 recombinant protein (Accession Number: XM\_047439733). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111519(HOMER3 antibody) at dilution of 1:500