

# Phospho-BCAR1-Y410 pAb



Catalog Number: 166519

## Product name

Phospho-BCAR1-Y410 pAb

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to Phospho-BCAR1-Y410

## Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding Y410 of human BCAR1

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - Phospho-BCAR1-Y410 pAb

Western blot analysis of extracts of HeLa cells, using Phospho-BCAR1-Y410 antibody at 1:5000 dilution. HeLa cells were treated by Pervanadate (1nM) for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 5s.