

Phospho-BCAR1-Y410 pAb



Catalog Number: 166519

Product name

Phospho-BCAR1-Y410 pAb

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to Phospho-BCAR1-Y410

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding Y410 of human BCAR1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-BCAR1-Y410 pAb

Western blot analysis of extracts of HeLa cells, using Phospho-BCAR1-Y410 antibody at 1:5000 dilution. HeLa cells were treated by Pervanadate (1nM) for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 5s.