

CPNE8 antibody

Catalog Number: 109519

Product name

CPNE8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CPNE8 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human lung tissue, human kidney tissue. Observed molecular weight by Western-blot: 63 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_153634). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109519(CPNE8 Antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung using Catalog No:109519(CPNE8 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung using Catalog No:109519(CPNE8 antibody) at dilution of 1:50 (under 40x lens)