

# Human HSPB11/C1orf41/IFT25/HSPC034 (His tag) recombinant protein



Catalog Number: 519519

## General Information

### Gene Name Synonym

Heat Shock Protein Beta-11, Hspb11, Placental Protein 25, PP25, HSPB11, C1orf41

### Protein Construction

A DNA sequence encoding the Human HSPB11/C1orf41/IFT25/HSPC034 Met1-Ser144 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 100mM NaCl, 2mM DTT, 10% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human HSPB11/C1orf41/IFT25/HSPC034 (His tag) recombinant protein