

# GDAP1 antibody



Catalog Number: 110919

---

## Product name

GDAP1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GDAP1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, fetal human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 36-41 kDa

## Preparation

This antibody was obtained by immunization of GDAP1 recombinant protein (Accession Number: NM\_001040875). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110919(GDAP1 antibody) at dilution of 1:500