

# MYCT1 antibody



Catalog Number: 112919

## Product name

MYCT1 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

MYCT1 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, BGC-823 cells, MKN-45 cells.

## Preparation

This antibody was obtained by immunization of MYCT1 recombinant protein (Accession Number: BC022556). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:112919(MYCT1 antibody) at dilution of 1:1000