

Anti-alpha Internexin Rabbit antibody



Catalog Number: 168919

Product name

Anti-alpha Internexin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

alpha Internexin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human alpha Internexin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M Nacl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

1:2000-1:10000

1:50-1:200

1:50-1:200

Validations



Western blot detection of alpha Internexin in HeLa cell lysates using alpha Internexin Rabbit mAb(1:1000 diluted).Predicted band size:55kDa.Observed band size:55kDa.