

EGR1 antibody



Catalog Number: 110219

Product name

EGR1 antibody

Specificity

human, mouse; other species not tested.

Antibody description

EGR1 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, Jurkat cells. Observed molecular weight by Western-blot: 68-70 kDa

Preparation

This antibody was obtained by immunization of EGR1 recombinant protein (Accession Number: NM_001964). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110219(EGR1 antibody) at dilution of 1:800