

# ABCB11 Polyclonal Antibody



Catalog Number: 162219

## Product name

ABCB11 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ABCB11

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1050-1150 of human ABCB11 (NP\_003733.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Western blot - ABCB11 Polyclonal Antibody


Western blot analysis of extracts of rat liver, using ABCB11 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 1s.

Immunofluorescence - ABCB11 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using ABCB11 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

# ABCB11 Polyclonal Antibody



Catalog Number: 162219

---

Immunofluorescence - ABCB11 Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using ABCB11 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.