

# Anti-OR5M3 antibody



Catalog Number: 176219

## Product name

Anti-OR5M3 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to OR5M3

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or5m3. aa range:53-102

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF

## Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

## Validations




Fig1.; Western Blot analysis of HEPG2 cells using Olfactory receptor 5M3 Polyclonal Antibody diluted at 1: 1000