

# Anti-PPP1A/PPP1CA+PPP1CB Rabbit antibody



Catalog Number: 167219

## Product name

Anti-PPP1A/PPP1CA+PPP1CB Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PPP1A/PPP1CA+PPP1CB Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human PPP1CA+PPP1CB

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/2000-1/10000


IHC: 1/20

ICC/IF: 1/50

FC: 1/20

IP: 1/20

## Validations



Western blot detection of PPP1CA+PPP1CB in K562,Rat Brain cell lysates using PPP1CA+PPP1CB Rabbit pAb(1:1000 diluted).Predicted band size:38kDa.Observed band size:38kDa.