

P-cadherin antibody

Catalog Number: 113619

Product name

P-cadherin antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

P-cadherin Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Positive WB detected in A431 cells, K-562 cells. Positive IHC detected in human breast cancer tissue. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 120 kDa

Validations




Preparation

This antibody was obtained by immunization of P-cadherin recombinant protein (Accession Number: NM_001793). Purification method: Antigen affinity purified.

A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113619(CDH3 antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000


Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113619(CDH3 antibody) at dilution of 1:100 (under 10x lens)



P-cadherin antibody


Catalog Number: 113619

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113619(CDH3 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-CDH3 (IP:Catalog No:113619,

4ug; Detection:Catalog No:113619 1:300) with A431 cells lysate 3040ug.



Immunofluorescent analysis of A431 cells using Catalog No:113619(CDH3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)