

UNC5D-Specific antibody



Catalog Number: 116619

Product name

UNC5D-Specific antibody

Specificity

human, mouse; other species not tested.

Antibody description

UNC5D-Specific Rabbit Polyclonal antibody.
Positive WB detected in mouse kidney tissue.
Observed molecular weight by Western-blot:
59-65kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_080872).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116619(UNC5D-Specific antibody) at dilution of 1:500