

CSTF3 antibody

Catalog Number: 109619

Product name

CSTF3 antibody

Specificity

Human; other species not tested.

Antibody description

CSTF3 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A431 cells. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 77-83 kDa

Preparation

This antibody was obtained by immunization of CSTF3 recombinant protein (Accession Number: NM_001033505). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109619(CSTF3 Antibody) at dilution of 1:600



IP Result of anti-CSTF3 (IP:Catalog No:109619, 4ug; Detection:Catalog No:109619 1:500) with A431 cells lysate 1200ug.