

GYG2 antibody



Catalog Number: 111319

Product name

GYG2 antibody

Specificity

human; other species not tested.

Antibody description

GYG2 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HEK-293 cells, human liver tissue. Observed molecular weight by Western-blot: 45 kDa, 66-70kDa

Preparation

This antibody was obtained by immunization of GYG2 recombinant protein (Accession Number: XM_011545600). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111319(GYG2 antibody) at dilution of 1:500