

Mouse ROBO4 (His Tag) recombinant protein



Catalog Number: 502319

General Information

Protein Construction

A DNA sequence encoding the mouse ROBO4 (NP_083059.2) N-terminal fragment (Met 1-Glu 480) was expressed, with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 38

Molecular Mass

The secreted recombinant mouse ROBO4 (aa 1-480) comprises 454 amino acids and has a calculated molecular mass of 49.6 kDa. As a result

of glycosylation, the apparent molecular mass of the recombinant protein is approximately 65 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse ROBO4 Protein (His Tag) SDS-PAGE