

CIRBP antibody

Catalog Number: 109319

Product name

CIRBP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CIRBP Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A549 cells, Cold treated HepG2 cells, HEK-293 cells, human brain tissue, mouse brain tissue, mouse testis tissue, PC-13 cells, Y79 cells. Positive IP detected in HepG2 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human breast cancer tissue, human placenta tissue. Observed molecular weight by Western-blot: 19 kDa

Preparation

This antibody was obtained by immunization of CIRBP recombinant protein (Accession Number: NM_001280). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IP, IF, WB

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations





WB result of CIRBP antibody (Catalog No:109319, 1:500) with si-control and si-CIRBP transfected



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109319(CIRBP antibody) at dilution of 1:1000



CIRBP antibody

Catalog Number: 109319



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109319(CIRBP Antibody) at dilution of 1:600

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109319(CIRBP antibody) at dilution of 1:100 (under 10x lens)



IP Result of anti-CIRBP (IP:Catalog No:109319, 3ug; Detection:Catalog No:109319 1:500) with HepG2 cells lysate 6000ug.

Immunofluorescent analysis of HepG2 cells using Catalog No:109319(CIRBP Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG