

FPR2 antibody

Catalog Number: 110719

Product name

FPR2 antibody

Immunogen

[Human FPR2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FPR2 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue.

Preparation

This antibody was obtained by immunization of FPR2 recombinant protein (Accession Number: NM_001462). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110719(FPR2 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110719(FPR2 Antibody) at dilution of 1:200 (under 40x lens).