

MSMB Polyclonal Antibody



Catalog Number: 161719

Product name

MSMB Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to MSMB

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-114 of human MSMB (NP_002434.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - MSMB Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MSMB antibody . Blue: DAPI for nuclear staining.