

STK32C antibody



Catalog Number: 115719

Product name

STK32C antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

STK32C Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, A549 cells, human skeletal muscle tissue. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of STK32C recombinant protein (Accession Number: BC015792). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115719(STK32C antibody) at dilution of 1:500