

Anti-NGFR antibody



Catalog Number: 176719

Product name

Anti-NGFR antibody

Specificity

Human, Mouse, Rat

Antibody description

Mouse monoclonal antibody to NGFR

Preparation

This antigen of this antibody was amino acids 393-427 at the c-terminus of ngfr p75 of human origin.

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, 40%Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IP, IF

Dilutions

WB: 1:100-1:1,000

IP: 1-2 µg per 100-500 µg of total protein(1 ml of cell lysate)

IF: 1:50-1:500

Validations




Fig1: Westernblot analysis of NGFR p75 expression in SK-N-NC whole cell lysate.




Fig2: Immunofluorescence staining of methanol-fixed HeLa cells showing membrane localization.