

# Anti-NGFR antibody



Catalog Number: 176719

## Product name

Anti-NGFR antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Mouse monoclonal antibody to NGFR

## Preparation

This antigen of this antibody was amino acids 393-427 at the c-terminus of ngfr p75 of human origin.

## Formulation

Liquid, 1\*TBS (pH7.4), 1%BSA, 40%Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IP, IF

## Dilutions

WB: 1:100-1:1,000

IP: 1-2 µg per 100-500 µg of total protein(1 ml of cell lysate)

IF: 1:50-1:500

## Validations




Fig1: Westernblot analysis of NGFR p75 expressionin SK-N-NC whole cell lysate.




Fig2: Immunofluorescence staining of methanol-fixed HeLa cells showing membrane localization.