

Anti-MSRB3 Rabbit antibody



Catalog Number: 168719

Product name

Anti-MSRB3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

MSRB3 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human MSRB3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

1:1000

Validations



Western blot detection of MSRB3 in C6 cell lysates using MSRB3 Rabbit mAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:21kDa.