

TGM2 antibody

Catalog Number: 116029

Product name

TGM2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TGM2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, rat heart tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 77 kDa, 82 kDa

Preparation

This antibody was obtained by immunization of TGM2 recombinant protein (Accession Number: NM_198951). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116029(TGM2 antibody) at dilution of 1:400



IP Result of anti-TGM2 (IP:Catalog No:116029, 4ug; Detection:Catalog No:116029 1:400) with HeLa cells lysate 1200ug.